

## SPEC SHEET

# Vertical Capsule Sorting Polisher

LSP-A/B MODELS

## INTRODUCTION

Capsule polishers and sorters are used in the pharmaceutical industry to clean capsules of their dust.

The dust collector is contained inside the body of the machine polisher meaning you do not need to connect an additional cleaning device.

**KEYWORDS** Capsule polisher | **Sorting** | Capsule cleaning | **Capsule dust** | **Auxiliary equipment** | Accessories





## ATTRIBUTES

- Vertical machine ensures better ergonomics and usage.
- Higher outlet meaning easier to connect to other devices such as a metal detector.
- Optional attachable vertical capsule sorting device to resolve capsule inconsistencies.

## TECHNICAL PARAMETERS

Model	LSP-A	LSP-B
Production Capacity	300,000 Capsules/Hour (1# Capsules)	600,000 Capsules/Hour (1# Capsules)
Air Consumption	22m <sup>3</sup> /h	44m <sup>3</sup> /h
Weight	100kg	
Dimensions	750×700×1450mm	
Power	0.65kw	
Power Supply	AC 380V 50Hz	
Air Pressure	≥0.3Mpa	



Whilst the information has been carefully compiled to the best of our knowledge, nothing is intended as representation or warranty on our part and no statement herein shall be construed as recommendation to infringe existing patents. In the effort to improve our products, we reserve the right to make changes judged to be necessary.